

# Sprout® PLUS Mini-Centrifuge

---

## Sprout® Plus Mini Centrifuge

With an ergonomic design and variety of colors, the Sprout® Plus Mini Centrifuge by Heathrow Scientific is loaded with features that bring increased efficiency to the lab

Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use

The lid was functionally designed to fit in the palm of the hand for easy repeat opening and quick tube extraction

A weighted rubberized base eliminates potential movement and dampens vibration and noise

The Sprout Plus is ideal for microfilter cell separations and HPLC samples, with adapters to accommodate all major microtube and PCR tube styles (1.5/2.0, 0.5, 0.2 mL) for customized use

The Sprout and Sprout Plus centrifuge housings are made from ABS plastic and lids are polycarbonate plastic enabling the unit to be easily cleaned

The Sprout Plus has the same footprint as the original Sprout, with a few added features and benefits:

A locking lid for added protection when the unit is in use

The lid can be unlocked and opened with one hand

A larger lid to accommodate longer and skirted tubes, as well as taller tube lids

A breaking function to quickly stop the unit when powered off or when the lid is opened

Meets the latest applicable North American and EU safety requirements

## Specifications:

Maximum Speed: Approx. 6,000 rpm (fully loaded)

Maximum RCF: 2000 Xg

Power: Suitable for worldwide voltages 100 to 240 VAC 50/60 HZ 0.6A Approx 7.3 W

Weight (w/round rotor): 1.60 lb, 0.73kg

Power: Input: 100-240VAC 50/60Hz, 1.0A, ) Output: 12VDC, 2.2A w/4 plug adapters for: UK, EU, North America, and Australia

For Indoor use, Double insulated, Level VI efficiency

Certifications: UL (Canada and US), GS, TUV, CE, PSE, and C-Tick

## Package includes:

Sprout Plus Mini Centrifuge  
6 place 1.5/2.0mL tube rotor  
16 place 0.2mL tube rotor (single or 2x8 strip)  
6 tube adapters 0.2mL  
6 tube adapters 0.5mL  
4 O-rings  
Storage case for rotors and adapters, includes a tube rack  
Low voltage, double insulated power adapter with 4 interchangeable plugs



