



Your **BROWNE** Process Indicator Autoclave Tape is now LEAD FREE

CURRENT





Clear result interpretation
 →
 New colour change
 Green to Brown/
 Black

NEW LEAD FREE



Clear result interpretation
 →
 New colour change
 Green to Brown/
 Black



Same roll size
 ■
 Same dimensions
 →
 New colour change
 ■
 Lead free

Current Order Information

CURRENT Steam Process Indicator Tape REF 0142 18mm (3/4") x 50m 48 rolls per carton	Steam Hi Tack Indicator Tape REF 0147 18mm (3/4") x 50m 48 rolls per carton
REF 0143 24mm (1") x 50m 36 rolls per carton	REF 0148 24mm (1") x 50m 36 rolls per carton

New Order Information

NEW Steam Process Indicator Tape REF 0162 18mm (3/4") x 50m 48 rolls per carton	Steam Hi Tack Indicator Tape REF 0167 18mm (3/4") x 50m 48 rolls per carton
REF 0163 24mm (1") x 50m 36 rolls per carton	REF 0168 24mm (1") x 50m 36 rolls per carton



Albert Browne Ltd, Chancery House,
 190 Waterside Road, Leicester. LE5 1QZ. UK
 Tel: +44(0)116 276 8636 ■ Fax: +44(0)116 276 8639
 Email: info@thebrownegroup.co.uk

■ The Company, whose policy is one of continuous improvement, reserves the right to make, without prior notification, changes in materials, dimensions and designs which, having regard to all circumstances, it thinks reasonable or desirable provided that the customer shall not be required to accept goods substantially different from those which are specified within the order.

■ NB: Colours shown are for reference only, and are as faithful as print processes allow.