




DURAN® HPLC bottle

with DIN thread, GL 45

With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete system comprising: DURAN® pressure plus laboratory bottle with 4-port screw cap (PP); four screw connections (black, M8 thread) and silicone seals. Connection of different hose diameters (1.6 mm and 3.2 mm) as well as sterile pressure equalisation sets (membrane filter 0.2 µm) is possible. Unused ports can be sealed with silicone blanking seals.

Typical application:

safe transfer of liquid media within a closed and sterile system (evaporation is reduced).



Ordering information:

Cat. No.	Capacity ml	DIN thread GL	d mm	h mm	Pack Quantity
11 298 21 *	500	45	86	176 **	2
11 298 20 *	1000	45	101	225 **	2

Accessories:

Cat. No.	Description	Pack Quantity
11 298 12	Screw cap HPLC, GL 45, 4 ports complete***	2
11 298 13	Spare part**** set for HPLC screw cap	1
11 378 01	Pressure equalization set 4-port cap (0.2 µm incl. Membrane filter)	1
11 298 19	Spare membrane filter for pressure equalization set, 0.2 µm	2

You can obtain these products from your specialist labware distributor.

* Replacement bottle.

** Bottle only.

*** Including silicone seals.

**** GL 45 screw cap, 4 x M8 screw cap, 12 x silicone seal.