

DURAN® MULTIPOINT CONNECTOR CAPS GL 32



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® MULTIPOINT CONNECTOR CAPS GL 32

2, 3 or 4 port connector caps for GL 32 DURAN® bottles

Applications

- Low volume sampling
- Media supply for small bioreactors
- Fed-batch bioreactor
- Perfusion circuit for cell culture
- Oxygenation / aeration of very small samples
- Use in a chemistry micro-reactor

Features

- Caps have a GL 32 thread that fit the smaller 50 mL round or 100 mL DURAN® square bottles
- Barbed male hose connector are suitable for tubing with a 3.2 mm (1/8") inner diameter (such as PharMed® BPT Tubing*)
- Available with 2, 3, or 4 surgical grade stainless steel (316L) barbed tubing ports
- Construction materials are certified as approved for food contact
- High temperature resistance (up to 180 °C)
- Autoclavable (@ 121 °C / 15 minutes)
- Reusable with a long service life

*PharMed® is a registered trademark of Saint-Gobain



Ordering information

	Description	Pack Unit	Cat. No.
NEW	DURAN® 2-port Connector caps GL 32, red PBT screw cap, PTFE insert and 2 stainless steel barbs	1	29 260 19 01
NEW	DURAN® 3-port Connector caps GL 32, red PBT screw cap, PTFE insert and 3 stainless steel barbs	1	29 261 19 02
NEW	DURAN® 4-port Connector caps GL 32, red PBT screw cap, PTFE insert and 4 stainless steel barbs	1	29 262 19 03

Recommended compatible DURAN® bottles

Description	Capacity (mL)	Pack Unit	Cat. No.
DURAN® Laboratory Bottle GL 32, clear, without screw cap. Retrace coded.	50	10	21 801 17 04
DURAN® Laboratory Bottle GL 32, amber, without screw cap. Retrace coded.	50	10	21 806 17 09
DURAN® Protect Laboratory Bottle GL 32, clear, safety coated, without screw cap. Retrace coded.	50	10	10 926 77
DURAN® Laboratory Bottle Square GL 32, clear, without screw cap. Retrace coded.	100	10	21 820 24 04

These products are available from your specialist laboratory dealer.



DWK Life Sciences GmbH
Hattenbergstrasse 10
55122 Mainz
Germany

Phone: +49 6131 - 1445 4131
Fax: +49 6131 - 1445 4016
sales@DWK-LifeSciences.com
www.DWK-LifeSciences.com

