



Pyrex® Centrifuge tube conical, ASTM

ASTM compliant 100ml conical centrifuge tubes

Manufactured to comply with a range of American standards and tests for water and sediment in petroleum products, these conical centrifuge tubes are manufactured from Pyrex® borosilicate glass.

- Complies with ASTM D96, D91, D-893 and D-1796
- In accordance with standard methods of testing for water and sediment in petroleum products
- Manufactured from chemically resistant Pyrex® borosilicate glass



Product Specification

Item Code	Approx overall length mm	Capacity ml	Neck ID mm	Pack quantity
3450/02MP	200	100	17	2