



DIN ISO 7550. DE-M marking for volumes greater than 5 µl. Calibrated to contain (TC, In). Color code to ISO for clear identification. CE-marked according to IVD-Directive 98/79 CE.

Minimum order quantity / quantity interval

box with 250 pieces

Variations of the product

Cat. No.	mark(s) at µl		color-code	length mm	A (%) ± %	reproducibility ≤ %
708707	1/2/3/4/5	*	white	125	0,30	0,6
708709	10		orange	125	0,25	0,5
708718	20		black	125	0,25	0,5
708722	25		2 x white	125	0,25	0,5
708727	40		2 x red	125	0,25	0,5
708728	20 + 40		2 x red	125	0,25	0,5
708733	50		green	125	0,25	0,5
708745	50 + 100		blue	125	0,25	0,5
708744	100		blue	125	0,25	0,5
708757	200		red	125	0,25	0,5