



**1621A Volumetric flasks with One graduation mark with glass/PE stopper ISO1042, DIN12664**

<b>Cat No</b>	<b>Cap. MI</b>	<b>Tolerance ±ml</b>	<b>Glass stopper</b>
1621A-01	5	0.025	7/11
1621A-02	10	0.025	7/11
1621A-03	25	0.040	10/13
1621A-04	50	0.060	10/13
1621A-05	100	0.100	12/14
1621A-06	200	0.150	14/15
1621A-07	250	0.150	14/15
1621A-08	500	0.250	16/16
1621A-09	1000	0.400	19/17
1621A-10	2000	0.600	24/20