


## Burettes, PTFE key, class AS, works certified

**PYREX®**

- Complies with ISO 385
- Individually tested and supplied with certificate showing exact volume dispensed at 3 points on the scale
- Individual serial number and certificate provides complete traceability to National Standards
- Each burette has a +30s waiting time before final reading can be taken
- Manufactured from precision-bore borosilicate glass tubing for high accuracy
- 1.5mm bore PTFE key which requires no lubrication
- White enamel graduations

Product Code	Capacity (ml)	Grads. ml	Tol. ± ml	Pack qty.
3295/22ASWC	10	0.02	0.02	1
3295/24ASWC	25	0.10	0.05	1
3295/26ASWC	50	0.10	0.05	1
3295/28ASWC	100	0.20	0.10	1



## Burettes, PTFE key, class AS

**PYREX®**

- Complies with ISO 385
- Manufactured from precision-bore tubing (borosilicate glass) for high accuracy
- Each burette has a +30s waiting time before final reading can be taken
- Graduations in white enamel for increased visibility with dark solutions
- Grease-free interchangeable PTFE key with 1.5mm bore

Product Code	Capacity (ml)	Grads. ml	Tol. ± ml	Pack qty.
3295/12AS	10	0.02	0.02	1
3295/14AS	25	0.10	0.05	1
3295/16AS	50	0.10	0.05	1
3295/18AS	100	0.20	0.10	1