

TEST TUBES / AMPOULES

CAPILLARY TUBES

Capillary Tubes

Soda glass

For melting point determinations.

Outside diameter 1.8/2.0mm, wall thickness 0.28/0.32mm.

Packed in tubular cartons of 100

Code	Length mm
G119/01	75
G119/02	100
G119/31	100 Sealed one end
G119/32	100 Sealed both ends



CENTRIFUGE TUBES

Centrifuge Tube

Soda glass

- Capacity 15ml

- Max. RCF x g: 1300

Code	Size mm	Wall thickness mm	Quantity per pack	Type
G131/01	112 x 17	1.0	100	Conical with rim
G131/71	112 x 17	1.0	100	Conical rimless
G132/01	112 x 17	1.0	100	Conical, graduated with rim
G138/01	110 x 17	1.0	100	Cylindrical with rim



FREEZE DRYING AMPOULES

Freeze Drying Ampoules

Neutral glass

Code	Capacity ml	Size mm	Quantity per pack
G225/01	0.5	100 x 7	625
G225/02	2.5	154 x 15	250

