



1601 Measuring Cylinder with spout & Round Base DIN 12680

Cat No	Capacity	Sub. Division(ml)	Tolerance <u>±</u> ml
1601-01	5	0.1	0.1
1601-02	10	0.2	0.2
1601-03	25	0.5	0.5
1601-04	50	1.0	1.0
1601-05	100	1.0	1.0
1601-06	250	2.0	2.0
1601-07	500	5.0	5.0
1601-08	1000	10.0	10.0
1601-09	2000	20.0	20.0