



Glass beakers, with spout, marking area and printed graduation

The Beaker is generally cylindrically shaped. It is flat bottom and with graduation. The material is Borosilicate glass 3.3. Most Beakers have spout which is easily for pouring. The below low form, tall form are widely used.

1101		Glass Beak low form, With spout and printed graduations	
Capacity(ml)	O.D.(mm)	Height(mm)	
5ml glass beaker	22	30	
10ml glass beaker	26	35	
25ml glass beaker	34	50	
50ml glass beaker	42	60	
100ml glass beaker	51	70	
150ml glass beaker	60	80	
200ml glass beaker	65	88	
250ml glass beaker	70	95	
300ml glass beaker	80	110	

400ml glass beaker	80	110
500ml glass beaker	87	118
600ml glass beaker	90	125
800ml glass beaker	100	135
1000ml glass beaker	106	145
2000ml glass beaker	130	185
3000ml glass beaker	150	210
5000ml glass beaker	170	270
10000ml glass beaker	217	350