



Cylinders, Nessler

MBL®

- For colour comparison of water samples etc.
- Manufactured from borosilicate tubing specially selected for freedom from colour or visible defects
- Cylinder bases are formed by a technique which gives maximum clarity
- Supplied in matched pairs for uniformity

Product Code	Capacity (ml)	Grads. ml	Tol. ± ml	O.D. x height mm	Pack qty.
CRC320	50	50	0.40	26 x 150	2
CRC322	100	50 & 100	0.80	34 x 180	2



Receivers, Crow, class B

PYREX®

- Complies with BS 658
- For the collection of distillate in oil and tar tests
- Graduations in white enamel
- Manufactured from Pyrex® borosilicate glass

Product Code	Capacity (ml)	Grads. ml	Overall height mm	Pack qty.
3400/04	100	1	250	2

LIPSOL® LABWARE DETERGENT

- Excellent cleaning performance - ideal for use with laboratory glassware and plasticware
- Great value for money
- Totally biodegradable

For full details see page 89

