

SPECIMEN TUBES SODA GLASS

Specimen Tubes

Flat bottom without closure

Code	Size mm	Wall thickness mm	Quantity per pack
G040/43	22 x 10	0.6	250

Specimen Tubes

Flat bottom supplied with polythene stoppers

Code	Size mm	Wall thickness mm	Quantity per pack
G050/17	38 x 10	0.6	250
G050/18	50 x 10	0.6	250
G050/19	38 x 12	0.85	100
G050/20	50 x 12	0.85	100
G050/22	75 x 12	0.85	100
G050/23	100 x 12	0.85	100
G050/24	50 x 16	0.85	100
G050/25	75 x 16	0.85	100
G050/26	50 x 19	0.85	100
G050/27	75 x 19	0.85	100
G050/28	50 x 25	0.85	100
G050/29	75 x 25	0.85	100
G050/30	100 x 25	0.85	100
G050/31	150 x 25	0.85	100

Spare Polythene Stoppers

Code	Size mm	Wall thickness mm	Quantity per pack
G051/00	10		1000
G051/01	12		1000
G051/02	16		1000
G051/03	19		1000
G051/04	25		1000



AGGLUTINATION TUBES

Dreyers Agglutination Tubes

Soda glass with funnel top

Code	Size mm	Wall thickness mm	Quantity per pack
G114/21	60 x 7	0.7	100

Dreyers Agglutination Tubes

Neutral resistance glass with funnel top

Code	Size mm	Wall thickness mm	Quantity per pack
G115/21	60 x 7	0.7	100

