


**KOCHGLÄSER UND ALLGEMEINES LABORGLAS**  
**BOILING FLASKS AND GENERAL LABORATORY GLASSWARE**

**DURAN® KRISTALLISIERSCHALE**  
**mit und ohne Ausguss**

*DURAN® CRYSTALLIZING DISH*  
*with and without spout*




Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Mit Ausguss DIN 12 338 With spout DIN 12 338				
21 311 24 <sup>1</sup>	20	40	25	10
21 311 32 <sup>1</sup>	40	50	30	10
21 311 34 <sup>1</sup>	60	60	35	10
21 311 38	100	70	40	10
21 311 41	150	80	45	10
21 311 44	300	95	55	10
21 311 49	500	115	65	10
21 311 54	900	140	75	10
21 311 59	2 000	190	90	10
21 311 63	3 500	230	100	10
Ohne Ausguss DIN 12 337 Without spout DIN 12 337				
21 313 24 <sup>1</sup>	20	40	25	10
21 313 32 <sup>1</sup>	40	50	30	10
21 313 34 <sup>1</sup>	60	60	35	10
21 313 38	100	70	40	10
21 313 41	150	80	45	10
21 313 44	300	95	55	10
21 313 49	500	115	65	10
21 313 54	900	140	75	10
21 313 59	2 000	190	90	10
21 313 63	3 500	230	100	10

<sup>1</sup> Ohne Aufdruck.

<sup>1</sup> Without print.

**DURAN® WÄGEGLAS**  
**mit eingeschliffenem Deckel**

**DURAN® WEIGHING BOTTLE**  
*with ground lid*



Durch den dichten Verschluss kann ein Probenverlust beim Transport nach dem Abwiegen verhindert werden. In niedriger und hoher Form erhältlich.

*Close-fitting lid prevents any sample loss during transport after weighing. Available in low and high forms.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d x h mm	Verpackungseinheit/ Stück Pack/Quantity
Niedrige Form <i>Flat shape</i>			
24 210 13	5	28 x 25	10
24 210 23	15	38 x 30	10
24 210 32	30	54 x 30	10
24 210 41	80	85 x 30	10
Hohe Form <i>Tall shape</i>			
24 211 13	10	28 x 40	10
24 211 18	20	32 x 50	10
24 211 23	45	38 x 70	10
24 211 24	70	44 x 80	10