



GHP-4

Uniform Heating at a Reasonable Price! This modern range of aluminium hotplate provides superior temperature uniformity. Chemical resistant top plate, easy to clean. Excellent heat conduction distributes the heat evenly and uniformly over the entire heating surface. Sensitive “stepless” temperature control is maintained by means of a hydraulic thermostat with minimum overshoot and undershoot Making possible the processing of delicate materials at low critical temperatures. Pilot lamp indicates when heat is on. Special radiation shields keep apparatus relatively cool during continuous operation. Supplied with a 3-wire line cord and sturdy screw type support device at the rear of the base to hold standard rod.

Temperature range: 40-300°C (optional 350°C). Top plate: 250x250mm.
Dimensions: W250xH150xD260mm. Watt: 1500W. Weight: 3.3kg.