

GMH-5

Giant Analog Aluminium Top Hotplate Magnetic Stirrer



Powerful Stirring and Uniform Heating! Quiet direct drive motor/magnet system. Chemical resistant top plate, easy to clean. Excellent heat conduction distributes the heat evenly & uniformity over the entire heating surface. Sensitive “stepless” temperature control is maintained by means of a hydraulic thermostat with minimum overshoot & undershoot Making possible the processing of delicate materials at low critical temperatures. Pilot lamp indicates when heat is on. Special radiation shields keep apparatus relatively cool during continuous operation. Supplied with a 3-wire line cord and sturdy screw type support device at the rear of the base to hold standard rod. Uniform, stable heating. Gentle or vigorous stirring.

Temp. Range: 40-300°C (optional 350°C). Speed range: 100-1400rpm. Max. capacity: 30liter. Top plate 250x250mm. Dimensions: W250xH150xD260mm. 1500W. Weight: 4.4kg.