

MH-1

Analog Aluminium Top Hotplate Magnetic Stirrer



Powerful Stirring and Uniform Heating! Quiet direct drive motor/magnet system. Chemical resistant top plate, easy to clean. Excellent heat conduction distributes the heat evenly & uniformity over the entire heating surface. Sensitive “stepless” temperature control is maintained by means of a hydraulic thermostat with minimum overshoot & undershoot. Making possible the processing of delicate materials at low critical temperatures. Pilot lamp indicates when heat is on. Special radiation shields keep apparatus relatively cool during continuous operation. Supplied with a 3-wire line cord and sturdy screw type support device at the rear of the base to hold standard rod. Uniform, stable heating. Gentle or vigorous stirring.