


Micro spoon spatulas



Material	A (Length)	a (Length)	B (Width)	D/Ø	F x E	Type	Article no.
18/10 Stainless	150	40	5	2,5	9x5	1	3340
18/10 Stainless	180	40	5	2,5	9x5	1	3341
18/10 Stainless	150	40	5	2	7x5	1	3342
18/10 Stainless	150	45	5	1,5	6x4	1	3343
18/10 Stainless	150	45	5	1,5	5x3	1	3344
18/10 Stainless	200	45	5	2,5	9x5	1	3345
18/10 Stainless	150	40	5	2	9x5	2	3350
18/10 Stainless	180	40	5	2	9x5	2	3351