

# SPATULAS AND SPOONS

# ספטולות וכפיות


- 2701- SPATULA - MICRO with solid handle, 23x2.5 mm blade.
- 2703- SPATULA POWDER - MICRO with solid handle, 23x3.2 mm blade.
- 2706- SPOON - MICRO with solid handle, 4.5 mm diam.
- 2710- WEIGHING SPATULA - SPOON 150,180 mm also in LEFT HAND.
- 2712- MICRO SPOON with SPATULA, Spoon: 7x10 mm.  
7x45 mm. In 150, 180, 210 mm lengths.
- 2715- SPOON - SPATULA 160, 180, 210 mm.
- 2720- WEIGHING SPOON, spoon size; 29 x 15 mm long.
- 2721- Same as 2720, but SEMI - MICRO, spoon size; 17x9 mm.  
130 mm long
- 2732- SPECIMEN SPATULA
- 2735- Same as 2732, but with wooden handle.
- 2742- SPATULA-NUFFIELD 150 x 10 mm. Made of 18/8 Stainless Steel.



2701

2703

2706



2710

2712

2715

2720

2721

2732

2735

2742