


## Micro spatulas, double end



Material	A (Length)	a (Length)	B (Width)	D/Ø	Article no.
18/10 Stainless	100	30	2	1	3010
18/10 Stainless	130	40	2	1	3011
18/10 Stainless	150	40	2	1	3012
18/10 Stainless	100	35	3	1,5	3013
18/10 Stainless	130	40	3	1,5	3014
18/10 Stainless	150	40	3	1,5	3015
18/10 Stainless	100	35	4	2	3016
18/10 Stainless	130	40	4	2	3017
18/10 Stainless	150	40	4	2	3018
18/10 Stainless	185	50	4	2	3019
18/10 Stainless	100	35	5	2	3020
18/10 Stainless	130	40	5	2	3021
18/10 Stainless	150	40	5	2	3022
18/10 Stainless	185	50	5	2	3023
18/10 Stainless	210	55	5	2	3024

18/10 Stainless	150	40	6	2	3025
18/10 Stainless	185	50	6	2	3026
18/10 Stainless	210	55	6	2	3027
18/10 Stainless	150	40	9	3	3028
18/10 Stainless	185	50	9	3	3029
18/10 Stainless	210	55	9	3	3030
Ni 99,5%	100	35	8	1,6	3046
Ni 99,5%	130	40	8	1,6	3047
Ni 99,5%	150	50	8	1,6	3053
Ni 99,5%	200	65	8	1,6	3057
PTFE coating	150	40	3	1,5	3710
PTFE coating	185	50	4	2	3711
PTFE coating	185	50	5	2	3712
PTFE coating	210	55	6	2	3713