



Flip & Grip™ Pipette Holder

The perfect choice for your lab's changing needs

This unique pipette holder gives you the option of using it on the bench top, or hanging it under a shelf to maximize available working space. Each Flip & Grip™ holder will accommodate up to 12 single channel pipettes, or a combo of 3 single and 3 multichannel pipettes. The unit holds all the most common brands of pipettes. Rubber grips prevent scratching of pipettes. The holder rotates for easy pipette selection and storage. Available in three colors.

- Weighted base. Fixing screw provided for under shelf option
- Measures 332 mm H x 193 mm W (bench top style); 187 mm L x 193 mm W (under shelf style)