

2134A**2132A**

Perfect for controlling sample temperature

These polypropylene racks keep microtubes afloat in a water bath or can be pushed into crushed ice. The detachable legs provide two options - attach legs so rack can stand upright on the bench, or detach legs so samples can be quickly popped out by pushing down onto any surface.

- Store down to temperatures of -80°C
- Material: Polypropylene
- Fits 1.5/2.0 mL tubes
- Color: White
- Sold in packs of 4



| Item No. | Well Shape | # of Wells | Well Diameter | Dimensions |
|----------|------------|------------|---------------|--|
| 2132A | Square | 16 | 10.9 mm | 8.4 x 4.2 x 2.6 in (21.34 x 10.67 x 6.6 cm) |
| 2134A | Round | 20 | 10.9 mm | 7.9 x 4.5 x 2.1 in (20.07 x 11.43 x 5.33 cm) |