



ONE-WELL CONNECTING TUBE RACKS 23052

Product Code:23052

- Made from polypropylene
- Interlock to cater for your personal requirements
- Form a line of 50 mL tubes, or use as a single, stand-alone rack.
- Molded circular base ring keeps conical or round-bottom tubes in place
- Tubes: 50 mL
- Wells: 1
- Well Dia: 29 mm
- 3.2 x 2.8 x 2.63 in (8.3 x 7 x 6.7 cm)
- 5/pk
- Assorted colors: blue, green, pink, yellow, orange
- Autoclavable



ONE-WELL TUBE RACKS (FIRCTION-FIT) 23053A

Product Code: 23053A

One-Well Tube Racks - Friction-Fit
Tube holders for small space requirements.

- Made from polystyrene
- Slides smoothly across the benchtop
- Tubes: 50 mL
- Wells: 1
- Well Dia: 29 mm
- 9.2 x 3 cm
- Not recommended for autoclaving