

Heavy Duty Vacuum Bottle

Heavy Duty Vacuum bottle

Material: PP Autoclavable

Code	Capacity ml
583255	2000
582255	4000



Filling/Venting Closure

Filling/Venting Closure

Autoclavable clouser for aseptic transfer of media, biological reagents, pure water and chemical to and from heavy duty vacuum bottle and carboys. Size 53 mm. can be used with Heavy duty vacuum bottle 2l ts. and 83 mm. can be used with 10 lts. and 20 lts. Carboy

Material: PP/silicone Autoclavable

Code	Capacity ml	Packing
583216	53	1
583217	83	1

