



Trough, multipurpose, polyethylene



- Rectangular robust "heavy-walled" trough
- Will stack inside each other for storage
- Capacity 3 litres

| Product Code | External dimensions mm | Pack qty. |
|--------------|------------------------|-----------|
| TGJ052 | 325 x 190 x 75 | 1 |



Trough, pneumatic, polyethylene



- Cylindrical, one piece pneumatic troughs
- Much lighter than glass
- Will withstand up to 80°C

| Product Code | O.D. x height mm | Pack qty. |
|--------------|------------------|-----------|
| TGJ058 | 200 x 100 | 1 |
| TGJ060 | 250 x 100 | 1 |
| TGJ062 | 300 x 150 | 1 |



Troughs, sterilising, polypropylene



- Rectangular sterilising troughs
- Robust and deep to accommodate a wide range of instruments etc.
- Can be autoclaved at 121°C

| Product Code | Capacity litres | l x w mm | Height without lid mm | Height with lid mm | Pack qty. |
|--------------|-----------------|-----------|-----------------------|--------------------|-----------|
| TWR202 | 1 | 210 x 122 | 78 | 91 | 1 |
| TWR204 | 10 | 390 x 290 | 110 | 130 | 1 |

Troughs lids, sterilising, polystyrene

- Lids are to be purchased separately and are not autoclavable

| Product Code | Description | Pack qty. |
|--------------|----------------|-----------|
| TWR203 | Lid for TWR202 | 1 |
| TWR205 | Lid for TWR204 | 1 |

