




### Aspirator, narrow neck, HDPE



- Large capacity aspirators with screw caps
- Ideal for bench-top dispensing of common solvents/distilled water etc. featuring smooth operating leak-proof stopcocks
- Good resistance to chemical attack

Product Code	Capacity litres	Neck I.D. mm	Overall height mm	Pack qty.
BGH034	5	49	330	1
BGH036	10	49	380	1
BGH038	25	80	580	1
BGH040	50	80	700	1

 Large carboys can be found on page 63.



### Aspirator, narrow neck, HDPE, spares


Product Code	Description	Pack qty.
SWQ016	Spare stopcock for BGH034 & BGH036	1
SWQ003	Spare stopcock for BGH038 & BGH040	1

### Aspirator, heavy duty, polypropylene



- Squat form for stable location on shelves, bench top or in transit
- Moulded handle for safe and secure carrying
- Clear moulded graduations allowing quick and easy volume assessment
- Leakproof heavy duty stopcock opens and closes in a quarter turn
- Autoclavable polypropylene for easy cleaning and sterilisation
- Caps have moulded drill guides for easy location of probes and feed-lines


Product Code	Capacity litres	Thread size mm	Max height x width mm	Pack qty.
BNP05A	5	45	270 x 220	1
BNP10A	10	90	360 x 280	1
BNP20A	20	90	420 x 350	1

 When autoclaving aspirators always ensure the caps are completely loosened or removed to prevent accidental implosion.

### Aspirators, heavy duty, polypropylene, spares

Product Code	Description	Pack qty.
CAP/045H	Spare cap for BNP05A/B	1
CGZ85PP	Spare cap for BNP10A/B & 20A/B	1
SCBNP01	Spare stopcock	1



 Carboys can be found on page 63.