



Beakers, straight sided, printed graduations, polypropylene



- Complies with ISO 7056
- Popular classic Azlon beakers
- Ideal for general laboratory use - straight stable shape
- Polypropylene material - good chemical resistance
- Can withstand temperatures up to 121°C and therefore autoclavable
- Graduations are clearly printed in indelible solvent resistant ink

Product Code	Capacity (ml)	Grads. ml	O.D. (mm)	Height mm	Pack qty.
BDA208P	10	1.0	32	32	10
BDA212P	25	2.0	42	40	10
BDA218P	50	5.0	51	51	10
BDA224P	100	10.0	62	62	10
BDA230P	250	25.0	82	82	10
BDA234P	400	50.0	94	95	10
BDA236P	500	50.0	93	101	10
BDA238P	600	50.0	106	107	5
BDA240P	1000	100.0	127	126	5
BDA242P	2000	200.0	158	160	1



i BDA234P and BDA238P are supplementary to ISO 7056.

Beakers, straight sided, printed graduations, polypropylene, sets



- Set of popular classic Azlon beakers
- Ideal for general laboratory use - straight stable shape
- Polypropylene material - good chemical resistance
- Can withstand temperatures up to 121°C and therefore autoclavable
- Graduations are clearly printed in indelible solvent resistant ink

Product Code	Description	Pack qty.
BDA218MES	Set of 3 beakers (50, 100 & 250ml)	1
BDA1198	Set of 3 beakers (500, 1000 & 2000ml)	1
BDA1199	Set of 5 beakers (10, 25, 50, 100 & 250ml)	1
BDA1200	Set of 8 beakers (10, 25, 50, 100, 250, 500, 1000 & 2000ml)	1
BDA300P	Set of 10 beakers (One of each size)	1