

FUNNELS, BUCHNER



Robust, lightweight. Unbreakable and resistant to corrosion and heat. Two piece construction. Top and bottom can be separated for easy and total cleaning. Autoclavable. Excellent chemical resistance. SUITABLE FOR FOODSTUFF.

Material:

PP

Art.	Filter Ø mm	Holes mm	Cap. ml	Stem length mm
437	42,5	1,2	40	47
438	55	1,1	70	57
439	70	2	180	57
440	80	2	285	65
441	90	2,5	390	68
442	110	2,5	810	92
443	160	2,75	2100	105
445	240	3,0	6000	143