

Desiccator stopcock

PYREX®

- For use with Pyrex® vacuum desiccators
- With Rotaflo® stopcock for vacuum tight seal


Product Code	Joint size	Pack qty.
1612/03	24/29	1



Desiccators, vacuum, polycarbonate

- Transparent polycarbonate cover and base
- Durable with high impact resistance
- Will safely withstand high vacuum pressures
- Ceramic plate permits the use of incandescent crucibles

Product Code	Description	I.D. mm	Plate size mm	Overall height mm	Effective inside height mm	Pack qty.
DWA215	With stopcock	215	205	262	178	1
DWA245	With stopcock	245	236	296	195	1
DWA246	Without stopcock	245	236	296	195	1



Desiccators, vacuum, polycarbonate, spares

Product Code	Description	Pack qty.
DWA205	Spare plate for DWA 215	1
DWA236	Spare plate for DWA245/246	1
DWA510	Spare stopcock for DWA215/245	1



Desiccator, cabinet, polystyrene





- Provides moisture-proof storage
- Clear high impact polystyrene cabinet
- Space saving stacking feature
- Airtight door with foam rubber gasket
- Supplied with sealed jar of silica gel

Product Code	Int. l x w x h mm	Ext. l x w x h mm	Door opening mm	Pack qty.	
DWA301	210 x 180 x 130	225 x 200 x 168	188 x 110	1	



Not suitable for vacuum use or filling with gases.

