

Hypodermic needles

Sterile, disposable needle in individual blister Luer Lock hub "Peel Pack". Ultra-thin wall and anti-friction treatment. Triple sharpened tip with "anti-coring" treatment International colour coding.

Indications for use:

Intradermal, subcutaneous, intramuscular and intra-articular injection; sampling of biological material, depending on the chosen gauge.

Colour	Gauge	MM.	Unit (pcs)
Yellow	20x1 ½	0,90x40	100
Black	22x1 ¼	0,70x30	100
Light blue	23x1 ¼	0,60x30	100
Light blue	23x1	0,60x25	100
Orange	25x5/8	0,50x16	100
Hazelnut	26x1/2	0,45x13	100
Grey	27x1/2	0,40x13	100
Green	21x1	0,80x25	100
Yellow	20x1	0,90x25	100
Cream	19x1	1,10x25	100
Green	21x1 ½	0,80x40	100
Black	22x1 ½	0,70x40	100
Cream	19x1 ½	1,10x40	100
Pink	18x1 ½	1,20x40	100

