


## CELLSTAR® Cell Culture Dishes



### Cell Culture Dishes

- Improved cell adhesion through physical surface treatment
- Vents ensure optimal gas exchange
- Sterile and user-friendly packaging
- Available in the sizes 35, 60, 94, 100 and 145 mm  $\phi$
- 8.7 to 143 cm<sup>2</sup> growth areas
- Easy stacking
- Maximal transparency for excellent microscopic analysis
- **New:** 35  $\phi$  dish also available with 4 internal wells

Like all Geiner Bio-One products, cell culture dishes are manufactured according to our high quality standards. Dishes are available in a wide variety of different dimensions and growth areas. In the case of 58 cm<sup>2</sup> and 143 cm<sup>2</sup> dishes, an extra high profile with a height of 20 mm is available.

For exact dimensions of our cell culture dishes, please refer to the product data sheets on our website.

Free of detectable DNase, RNase, human DNA  
non-pyrogenic

non-cytotoxic

Cat.-No.	627 160	627 170	628 160	632 171
$\phi$ [mm] x height [mm]	35 x 10	35 x 10	60 x 15	94 x 16
Growth area [cm <sup>2</sup> ]	8.7	0.93/well	21	58
Total volume [ml]	10	10	28	80
Working volume [ml]	5	0.1/well	6 – 7	15 – 16
Vents	+	+	+	-
TC surface treatment	+	+	+	+
Sterile	+	+	+	+
Quantity per bag/case	10/740	10/740	10/600	10/480

New

Cat.-No.	664 160	639 160
$\phi$ [mm] x height [mm]	100 x 20	145 x 20
Growth area [cm <sup>2</sup> ]	58	143
Total volume [ml]	100	240
Working volume [ml]	16 – 17	25 – 27
Vents	+	+
TC surface treatment	+	+
Sterile	+	+
Quantity per bag/case	15/360	5/120