



Mega Rack

High capacity storage for your tubes

Fully assembled, polypropylene Mega racks are great in clinical labs for storage of high quantity samples. Sturdy, light weight material for easy transport. The rack has identifying alpha-numeric for quick location and retrieval of tubes. Conveniently place different test groupings into one rack and eliminate need for multiple racks. Available in two sizes to fit standard 12/13 mm or 16 mm tubes. Choose from blue or white. Measures 337 x 267 x 91 mm. Eac

Cat No.	Description	Capacity	Color
HD24322A	10-13 mm Tube Rack	216 Tubes	Blue
HD24322E	10-13 mm Tube Rack	216 Tubes	White
HD24323A	13-16 mm Tube Rack	120 Tubes	Blue
HD24323E	13-16 mm Tube Rack	120 Tubes	White