

A**A - SCOOP WITH FLAT BOTTOM**

Material: Polypropylene
 Color: Neutral - For food use
 For withdrawal of solids and granulate samples.
 The measure includes with handle.

B - LARGE CUP

Material: High density Polyethylene
 Color: White - For food use

B

Item	Capacity ml	Total length mm	Useful length mm
PL 66	2000	400	305

Item	Capacity ml	Total length mm	Useful length mm
PL 2147	3000	200	115

SCOOPS

Material: High density Polyethylene
 Color: Neutral - For food use
 For withdrawal of solid and granulates samples.
 They have the rounded bottom for a better use.



Item	Capacity ml	Total length mm	Useful length mm
PL 62	90	180	125
PL 63	250	240	155
PL 64	600	300	195
PL 65	1100	360	240
PL 67	2500	400	285