

Features:

Uses polystyrene raw materials that meet USP-Class VI standards and is produced in a GMP

100,000 clean workshop

7 kinds of capacity specifications, 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Two packaging methods, independent paper-plastic packaging and bulk packaging

The end of the tube is marked with a color ring, and the body of the tube has a corresponding color bar background, which is convenient for quickly identifying the capacity

Positive and negative double scale design, clear printing, improve sampling accuracy

Negative scale design provides a certain additional sampling capacity

No pyrogen, no endotoxin, no cytotoxicity

Irradiation sterilization, SAL 10-6

Specifications:

Catalog number	capacity (ml)	color ring	Color code	package style	Sterilize	Branch/Pack	Bag/box	Box/carton	Boxing Quantity
311000	1		YES	individual package	+	100	5	2	1000
312000	2		YES	individual package	+	100	4	2	800
313000	5		YES	individual package	+	50	6	2	600
314100	10		YES	individual package	+	50	5	2	500
315100	25		YES	individual package	+	50	4	2	400
316100	50		YES	individual package	+	25	4	2	200
317110	100	-	NO	individual package	+	1	50	2	100
313100	5		YES	Welding type individual package	+	50	6	2	600
313211	5		NO	In bulk	+	25	20	6	3000
314311	10		NO	In bulk	+	25	14	6	2100