

Thermotech™ Beakers

- ✓ Isostatically molded in pure, FDA approved, PTFE
- ✓ Machined base for good heat transfer
- ✓ Molded-in pouring spout

Totally chemically resistant and inert



Beaker and Beaker Cover

Beaker Reference	Beaker Volume ml (Fl oz)	Beaker Height	Beaker Diameter	Cover Reference	Cover Diameter
007.001	1 (0.03)	18 mm	13 mm	013.020	20 mm
007.005	5 (0.17)	26 mm	20 mm	013.030	30 mm
007.010	10 (0.33)	33 mm	24 mm	013.030	30 mm
007.025	25 (0.83)	47 mm	32 mm	013.040	40 mm
007.050	50 (1.67)	60 mm	43 mm	013.050	50 mm
007.100	100 (3.33)	68 mm	54 mm	013.065	65 mm
007.150	150 (5.00)	69 mm	59 mm	013.075	75 mm
007.250	250 (8.33)	97 mm	66 mm	013.080	80 mm
007.400	400 (13.3)	106 mm	80 mm	013.100	100 mm
007.500	500 (16.7)	125 mm	80 mm	013.100	100 mm
007.600	600 (20.0)	125 mm	90 mm	013.125	125 mm
007.1000	1000 (33.3)	155 mm	100 mm	013.125	125 mm
007.2000	2000 (66.6)	210 mm	120 mm	013.150	150 mm
007.5000	5000 (167)	280 mm	170 mm	N/A	N/A

Thermotech™ PTFE Beakers



This product features the unique combination of a pure PTFE body with a specially formulated stabilised PTFE Carbon base as a single molding to give:

- ✔ Totally inert pure PTFE inner surface
- ✔ Heatable to 270°C without distortion
- ✔ Improved heat transfer

A Cowie Registered Design

Beaker Reference	Beaker Volume ml (Fl oz)	Beaker Height (mm)	Beaker Diameter (mm)
007.0100	100 (3.33)	74	56
007.0250	250 (8.33)	94	75
007.0400	400 (13.3)	112	85

