

cowie® PTFE Stirring Magnets

[HOME](#)
[LABWARE INDEX](#)


Hub

Very stable stirrer, especially at low speeds



Do-Nut

Designed specifically for automated systems. For use in tubes and allows the use of robotic or similar sampling directly through the stirrer magnet. Rare Earth magnet core.



Double-Ended

Double paddle action for efficient stirring plus high stability. Colour coded.

Note: coloured PTFE coatings are not as inert as pure PTFE

Reference	L (mm) x D (mm)	Reference	D (mm) x Ht x Hole (mm)	Reference	L (mm) x D (mm)
001.2301	45 x 27	001.3110	9.5 x 13 x 2	001.1335	35 x 8
001.2302	62 x 37	001.3113	12.5 x 15 x 3	001.1355	55 x 8
001.2303	70 x 37			For coloured use suffix: -R= Red -B = Blue -Y = Yellow e.g. 001.1335-Y	

NB: The dimensions of stirrer bars are nominal and approximate to the stated value by $\pm 5\%$ for length and $\pm 10\%$ for diameter