

HOME

LABWARE INDEX


cowie® Stirring Magnets



For profile drawing see:
Octohedral
Micro

Colour Coded

For use where identification is important. NOTE: Coloured PTFE coatings are not as inert as pure PTFE



Glass Covered

For use with very abrasive media which may erode PTFE

Reference	L (mm) x D (mm)	Reference	L (mm) x D (mm)
Octohedral		001.1206	6 x 5
001.513	Use suffix R = Red B = Blue Y = Yellow	001.1212	12 x 5
001.515		001.1225	25 x 6
001.522	e.g. 001.513-Y	001.1245	45 x 8
001.525		001.1260	60 x 8
001.538	See 001.500 Series for details		
001.551			
001.575			
001.1140			
001.1160			
Micro			
001.802	Use suffix R = Red B = Blue Y = Yellow		
001.803			
001.805	e.g. 001.802-R		
001.806			
001.807	See 001.800 Series for details		
001.808.3			
001.810			
001.813			
001.815			

NB: The dimensions of stirrer bars are nominal and approximate to the stated value by ±5% for length and ±10% for diameter

cowie®