

cowie[®] PTFE Tissue Grinders

[HOME](#)
[LABWARE
INDEX](#)

- ✓ PTFE, stainless steel and precision borosilicate glass
- ✓ Ultra pure - no contamination
- ✓ Standard clearance 0.15 mm - 0.25 mm

A wide variety of clearances available to order



Tissue Grinder / Homogeniser

Glass Vessel

Reference	Vol. ml (Fl oz)	Bore Dia.	Height
011.102	2	8 mm	120 mm
011.105	5	12 mm	135 mm
011.110	10	15 mm	150 mm
011.115	15	19 mm	155 mm
011.130	30	25 mm	175 mm
011.150	50	32 mm	195 mm

Volumes stated are working volumes with the plunger in place

Plunger - Plain

Plunger - Serrated

Reference	Vol. ml (Fl oz)	Height	Reference	Vol. ml (Fl oz)	Height
011.202	2	230 mm	011.302	2	230 mm

011.205	5	235 mm	011.305	5	235 mm
011.210	10	270 mm	011.310	10	270 mm
011.215	15	270 mm	011.315	15	270 mm
011.230	30	270 mm	011.330	30	270 mm
011.250	50	270 mm	011.350	50	270 mm

Volumes stated are working volumes with the plundger in place

cowie[®]