

[HOME](#)**cowie**<sup>®</sup> PTFE Thin Sleeves[LABWARE  
INDEX](#)

- ✓ Self-release seal to prevent sticking
- ✓ Operates at medium to low vacuum
- ✓ Use up to 280°C (540°F)

Eliminates the need for grease

**Thin Sleeves**

<b>'A' Length Cone - USA</b>		<b>'B' Length Cone - Europe</b>	
<b>Reference</b>	<b>Joint</b>	<b>Reference</b>	<b>Cone</b>
019.007	7 / 25	019.107	7 / 16
019.010	10 / 30	019.110	10 / 19
019.012	12 / 30	019.112	12 / 21
019.014.1	14 / 20	019.114	14 / 23
019.014	14 / 35	019.119	19 / 26
019.019.1	19 / 22	019.124	24 / 29
019.019	19 / 38	019.129	29 / 32
019.024	24 / 40	019.134	34 / 35
019.029	29 / 42	019.140	40 / 38
019.034	34 / 45	019.145	45 / 40
019.040	40 / 50	019.150	50 / 42
019.045	45 / 50	019.155	55 / 44
019.050	50 / 50	019.160	60 / 46
019.055	55 / 50	019.170	71 / 51
019.060	60 / 50	019.185	85 / 55
019.070	71 / 60		

**cowie**<sup>®</sup>